# Aprilla

ABS 1055

**Body Analysis Digital Bathroom Scale**

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
ES

PT

- *This device is designed and manufactured for home use. It is not suitable for commercial or industrial purposes.*
- *This device may be interest of children, but as you know it is not a toy. In this reason, device being where children can't reach place.*
- *Other information in this guide to see the section "**Safety Instructions**".*
- *Please read the manual, you can start using the product after that.*

#### ***Children Being Supervised not to Play with the Appliance***

*this appliance can be used by children aged from 8 years and above and persons with reduced physical, sensory or mental capabilities or lack of experience and knowledge if they have been given supervision or instruction concerning use of the appliance in a safe way and understand the hazards involved.*

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**EN** Dear Aprilla User,

*First of all thank you for choosing Aprilla product. Please read this manual carefully before using your device. Please keep the manual as it may need in future. This guide is designed to help you to achieve maximum performance.*

*We wish you a pleasant use.*

**Aprilla**

## SAFETY INSTRUCTIONS

EN

• Please read the instruction manual carefully before using the device. Please keep these instructions, the guarantee certificate, the sales receipt and, if possible, the carton with the inner packaging. • This appliance is not intended to be used by persons (including children) with reduced physical, sensory or mental capabilities, or lack of experience and knowledge, unless they have been given supervision or instructions concerning the use of the appliance by a person that is responsible for their safety. • By ignoring the safety instructions the manufacturer can not be held responsible for the damage. • To protect children against the dangers of electrical appliances, please make sure that you never leave the appliance unattended. Therefore you have to select a storage place for the appliance where children are not able to grab it. Make sure that the cable is not hanging in a downward position. • This appliance is only to be used for household purposes and only for the purpose it is made for. • The appliance must be placed on a stable, level surface. • Handle the balance with care, as it is a precise measuring instrument. Do not drop it or jump on it. • Do not load the appliance with more than 180kg to avoid damaging the device. • If the scales are carried or not going to be used, please store them only in a horizontal position. The body fat analysis scale is not suitable for the following persons: • Children under 10 years and adults over 85 years. • Persons with symptoms of fever, oedema, or osteoporosis. • Persons undergoing treatment for dialysis. • Persons who use a pacemaker. • Persons who take cardiovascular medication. • Pregnant women. • Sportspersons who perform over 5 hours of intensive sporting activities each day. • Persons with a heartbeat lower than 60 beats/minute. Batteries • Do not expose the battery to high temperatures or direct sunlight. Never throw batteries into the fire. There is a danger of explosion! • Keep batteries away from children. They are not a toy! • Do not open the batteries by force. • Avoid contact with metallic objects. (Rings, nails, screws et cetera.) there is a danger of short-circuiting! • As a result of a short-circuit batteries may heat up considerably or even catch fire. this may result in burns. • For your safety the battery poles should be covered with adhesive strips during transport. • If a battery leaks do not rub the liquid into contact with the liquid they should be washed thoroughly and eyes rinsed out with clear water, if symptoms persist please consult a doctor."

## BEFORE THE FIRST USE

• Take the appliance and accessories out the box. Remove the stickers, protective oil or plastic from the device. • The device can only function correctly if you stand on the contact surface with bare, dry feet! • The main principle of this device is based on the measurement of the electric impedance within the human body. In order to determine the actual body data, the device leads a weak, unnoticeable electric signal through your body. • This signal is generally safe! However, please read our special safety instructions! • The method of measurement is referred to as "bioelectric impedance analysis" (BIA), a factor which is based on the relationship between body weight and associated body water percentages as well as on other biological data (age, gender, body weight). • There are heavy persons with a well-balanced percentage between fat and muscle mass, and there are seemingly slim persons with

significantly high body fat values, fat is a factor that puts serious strain on the body, in particular the cardio-vascular system. Therefore, controlling both values (fat and weight) is a prerequisite for health and fitness.

## USE

• The scales are factory set with the units "kg" and "cm". To change the settings to "lb" or "kg" and "Ft:in", proceed as follows: On the back side of the scale beside the battery compartment a switch can be found, push the below switch, choose between the measuring units kg/lb/st.

## WEIGHING FUNCTION

• Place the scale on a hard, flat surface (avoid carpets). An uneven floor contributes to inaccurate measurements. • Stand on the scales with both feet. The device switches on automatically. • Distribute your weight equally and wait for the scale to determine your weight. • The digits on the display will flash before your exact weight will be shown. • In order to extend the battery life, the device will automatically turn off 10 seconds after you have stepped off the scale.

## SAVING PERSONAL DATE

• The scale has default values that you can change, as described below. • Carry out the following steps in sequence and within a few seconds, otherwise the device will automatically cancel the input mode. • Press the button SET. The memory space will flash on the display (P0). Confirm the setting by pressing SET or select another memory space by pressing the UP/DOWN button and then SET to go to the next setting. • The symbol for "male" or "female" will flash. Confirm the setting by pressing SET or use the UP/DOWN button to choose the symbol for your gender and then press to go the next setting. • The height setting will flash. Select your height by pressing the UP or DOWN button. Holding down one of the buttons will accelerate the process. Confirm your entry by pressing SET. • The setting AGE will flash. Select your age by pressing the UP or DOWN button. Holding down one of the buttons will accelerate the process. Confirm your entry by pressing SET. • The device can only function correctly if you stand on the device with bare, dry feet. If necessary wipe the contact surface to improve the results. • Due to the significant fluctuation of body fat the optimal testing time is between 6 p.m. and 8 p.m., since the body fat remains very consistent during this period of time. • Press the SET button, the user number will flash (P0-P9). Use the UP/ DOWN button to select your pretest data. • Wait for the user number to stop flashing. • Stand on the contact surfaces with both feet. • Firstly your weight will flash on the display. • After a few seconds, the following results will be shown: • the body fat content in% (FAT symbol) • the results analysis in the form of a bar • the body water content in % (TBW symbol) • the muscle mass content in % (MUS symbol) • the bone mass content in % (BON symbol)

## ERROR

If an error message appears on the display, wait until the device is turned off. Turn it on again to restart the operation. If the error display has not disappeared, take the batteries out of the battery compartment for 10 seconds. • 8888 = Self test (when

# Aprilla


switching on) to check if the numeric display functions properly. • `---- = The scale is not placed properly. Err = The percentage of body fat content lies above or below the limit of the measuring scale • EEEE = The weight lies above the allowed maximum weight.

## CLEANING AND MAINTENANCE

EN

• Do not immerse the appliance in water or any other liquid. • Clean the appliance with a damp cloth. Never use harsh and abrasive cleaners, scouring pad or steel wool, which damages the device.

DISPLAY SYMBOL	MEANING
Kg/lb/st)	Body weight
P0-9	User memory space
♂	Male
♀	Female
CM(FT;In)	Body height
Age	Age
FAT	Body fat display
TBW	Body water display
MUS	Muscle mass display
BONE	Bone weight display
%	Percentage-body fat/body water/muscle mass/bone weight

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**Do not operate device without reading the manual.**

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**This device appropriate for regulations according to Waste Electrical and Electronic Equipment (WEEE).**

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**This product is produced in nature friendly factory.**

© 2015 BMVA. All rights reserved. Aprilla logo is owned by BMVA and registered. All other trademarks and signs are property of their owners. BMVA is not responsible for printing and composing mistakes.